



DKS600

PRESIDENT

OWNER'S MANUAL

BLADEZ HEALTH PRODUCTS CORPORATION
20155 ELLIPSE
FOOTHILL RANCH, CA 92610

TELEPHONE: 1-866-610-0330
EMAIL: BLADEZHEALTHPRODUCTS@EBLADEZ.COM
WEBSITE: WWW.EBLADEZ.COM




TABLE OF CONTENT

1.	PREFACE.....	3
2.	ASSEMBLY INSTRUCTIONS.....	4
3.	SAFETY NOTICE.....	6
	3.1 Before Driving	
	3.2 While Driving	
4	ELECTROMAGNETIC INTERFERENCE.....	9
5.	PARTS INTRODUCTION.....	11
	Parts Description	
6.	OPERATION.....	12
	6.1 Operation Panel	
	6.2 How to Operate Your vehicle	
	6.3 How to Set to Freewheel Mode	
7.	CHARGING AND BATTERY.....	19
	7.1 Charging the Battery	
	7.2 Charger	
	7.3 Battery	
8.	INSPECTION AND MAINTENANCE.....	22
	8.1 Daily Checking	
	8.2 Regular Checking Record	
	8.3 Fuse and Tire	
	8.4 Maintenance	
	8.5 Storage	
9.	TROUBLE SHOOTING.....	26
10.	SPECIFICATIONS.....	28
11.	WARRANTY.....	29
	11.1 Vehicle Serial Number	
	11.2 Warranty Policy	

1. PREFACE

Please carefully read this owner's manual before using the vehicle. Improper use of the vehicle could result in harm, injury or accidents.

- ◆ This owner's manual includes instructions for all aspects of the vehicle, including assembly instructions and how to deal with possible accidents.
- ◆ This owner's manual is written solely for the Bladez **B*Mobile, DKS600 - President**
- ◆ The symbols used in this manual are explained below.

 Warning	Improper use of this vehicle could lead to serious injury
 Attention	Improper use of this vehicle could lead to personal injury and/or property damage.
 Suggestion	Follow these instructions to maintain safe and efficient operations of the vehicle.

- ◆ This manual includes a repair and maintenance record chart and the product registration card. Please keep them in a safe place.
- ◆ If someone else uses the vehicle, make sure that you give them the instruction manual for his or her reference.
- ◆ We reserve the right to make design modifications without notifying the buyer. It is your responsibility to maintain an updated manual either by requesting it through mail or visiting our website. (www.ebladez.com)

“Caution: Federal law restricts this device to sale by or on the order of a practitioner licensed by the law of the State in which he/she practices.”

2. ASSEMBLY INSTRUCTIONS

1. Once you have removed your vehicle from the box it will be folded. (Fig.1) The basket and the seat are included with it. Remove the seat and basket and set them aside.
2. Once you have removed the seat and basket, locate the tiller handle and fold the handlebars up into the desired position. (Fig. 2)

FIG. 1



FIG. 2



FIG. 4



FIG. 3



3. Before you install your seat you will need to connect the batteries. Remove the battery cover by removing the screw with a Phillips screw driver to gain access to the batteries. (Fig. 3) Connect the wires from the batteries to the wire coming from the controller. (Fig. 4) Once you have connected the batteries replace the cover.

4. Locate the basket, water bottle holder, and cell phone basket. Attach your water bottle holder and your phone basket to the large basket using the screws provided. (Fig 5 and 6) Once you have your basket ready you can now attach it to the vehicle.

FIG. 5




FIG. 6



2. ASSEMBLY INSTRUCTIONS

5. To attach the basket to the vehicle use the five screws provided. Place the basket on the metal bracket on the front of the vehicle and line up the five holes in the bottom of the basket (Fig. 7) to the holes on the metal bracket. Insert the screws and tighten to secure the basket.



6. Lastly, install the seat. Line up the pin on the bottom of the seat with the seat post and slide it into place. (Fig. 8) Lift the lever on the side of the seat to lock the seat into place. This lever will also be used to turn the seat in a 360 degree circle when needed.



3. SAFETY NOTICE

3.1 Before Driving

The user should become familiar with all functions and features of this vehicle before driving. Please always keep safety in mind. It may be beneficial to have an assistant such as a family member or a care professional present while using the vehicle. To minimize any personal injury please read all of the following safety precautions prior to driving.

➤ **The same traffic rules that apply to pedestrians apply to the use of this vehicle**

For your safety, please follow all rules that apply to pedestrians.

- ◆ Ride on the sidewalk, parkway or pedestrian areas only. Never ride on Freeways or roads. It may be difficult for traffic to see you when driving your vehicle.
- ◆ Always be aware of all obstacles and conditions around you. (Cars, people, terrain, and weather).
- ◆ Take into consideration the weather conditions before driving.
- ◆ Be extremely cautious when driving your vehicle in busy areas, such as shopping malls, business centers, etc.
- ◆ If you take prescription drugs please consult your physician before operating your vehicle. Some medications may impair your ability to maintain safe operations of your vehicle.
- ◆ Please do not operate while under the influence of alcohol, this may impair your ability to operate your vehicle in a safe manner.
- ◆ Please do not drive your vehicle at night.
- ◆ Please obey pedestrian traffic lights.

➤ **Practice operating your vehicle**

The user must use care and common sense when driving your vehicle. While becoming familiar with the operation of your vehicle, practice in a wide open area, like a park. In order not to fall out of your vehicle while driving, keep in mind the status of the driving motion, such as accelerating, stopping, turning, reversing, and up-and-down ramps.

- ◆ Please set the speed dial to the low setting during initial practice.
- ◆ Be sure you are able to control and operate your vehicle easily and confidently. Then you can adjust your vehicle to a higher speed

3. SAFETY NOTICE

➤ **The B*Mobile vehicle is limited to one passenger**

Do not carry passengers on your vehicle. (Including children)

➤ **Do not use this vehicle to carry or haul goods**

◆ The maximum weight can be carried (including goods) is 308.6 lbs (140 kg). Refer to section “10-SPECIFICATIONS” on page 28 “WEIGHT CAPACITY.”

◆ Maximum loading weight for the basket is 11 lbs (5kg).

3.2 While Driving

➤ **Do not stretch your body out of the vehicle while in motion**

◆ Such action may cause you to lose balance or to be injured due to falling.

◆ Pay attention to avoid your clothes being tangled into the wheels.

➤ **When Crossing The Train Tracks**

◆ Before entering the crossing please stop completely. While crossing please watch that the vehicle’s wheels do not get stuck in the track.

➤ **Do not use your vehicle under circumstances below.**

◆ On a road with heavy traffic, mud, gravel, bumps, snow or ice. Use caution on roads that have a canal or waterway without any fences, or hedges or that are too narrow. These conditions can affect the performance of your vehicle and could lead to an accident or injury.

◆ Do not drive at night or in rain, snow, or mist. Driving in these conditions may affect your ability to drive in a safe manner and could lead to an accident or injury

◆ Do not drive your vehicle in a zigzag motion or make any sudden turns. Sudden turns or driving in a zigzag motion can result in the tipping over of your vehicle and could lead to possible injury or damage to your vehicle. To avoid accidents, practice making turns to become familiar with how your vehicle maneuvers.

◆ Your vehicle is not designed to travel up or down stairs or escalators, always use an elevator.

➤ **Mobile phones and other electronic equipment**

◆ Do not use a mobile phone or other wireless communication device while driving.

◆ Do not charge the mobile phone or other electric device from your vehicle battery

3. SAFETY NOTICE

➤ **Ramps, inclines and declines**

- ◆ When coming to an inclined road, Set the speed dial to a middle setting and proceed slowly. Always maintain a safe speed and avoid any sudden starts or stops. When going down a declined road, set the speed dial to a low setting to ensure a safe and controlled decent. Never put the vehicle in reverse when traveling down a decline.
- ◆ Do not drive over large drop offs. For more details you can refer to section “10-SPECIFICATIONS” on page 28 “MAX GROUND CLEARANCE.”
- ◆ Do not cross wide gutters or you may get your wheels stuck.
- ◆ Please slow down when driving on roads with drops offs.
- ◆ Do not make sudden turns when driving on gravel roads or ramps.
- ◆ Avoid any terrain or obstacles that make you feel uncomfortable or if have any doubts that you vehicle will not perform in a safe manner.
- ◆ Many buildings now have ramps with 180 degree turns, so it is a good idea to practice your turning skill in order to not get stuck on railings or to get caught on any obstacles witch could lead to injury.

➤ **Maximum User Weight Limit**

Refer to section “10-SPECIFICATIONS” on page 28 “WEIGHT CAPACITY.” Loading over the limitation will lead to the damage of your vehicle or malfunction, which could be unsafe. The warranty does not cover this type of over-loading damage.

4. ELECTROMAGNETIC INTERFERENCE

This portion of the owner's manual will provide the user with basic information about the problems with EMI, known sources of EMI, protective measures either to lessen the possibility or exposure, or to minimize the degree of exposure. This section also shows some conditions when unexpected or erratic movements may occur.

Caution: It is very important that you read this information regarding the possible effects of electromagnetic interference on your electric B*Mobile vehicle.

➤ **ELECTROMAGNETIC INTERFERENCE (EMI) FROM RADIO WAVE SOURCES**

Powered vehicles may be susceptible to electromagnetic interference (EMI), which is interfering electromagnetic energy (EM) emitted from sources such as radio stations, TV stations, amateur radio (HAM) transmitters, two-way radios, and cellular phones. The interference (from radio wave sources) can cause the electric powered vehicle to release its brakes, move by itself, or move in unintended directions. It can also permanently damage the electric powered vehicle's control system. The intensity of the interfering EM energy can be measured in volts per meter (V/m). Each electric powered vehicle can resist EMI up to certain intensity. This is called its "immunity level". The higher the immunity level is, the greater the protection is. Current technology is capable of achieving at least a 20 V/m immunity level, which would provide useful protection from the more common sources of radiated EMI. This electric powered vehicle model as shipped, with no further modification, has an immunity level of 20 V/m without any accessories.

There are a number of sources of relatively intense electromagnetic fields in our everyday environment. Some of these sources are obvious and easy to avoid. Others are not apparent and exposure is unavoidable. However, we believe that by following the warning listed below, your risk to EMI will be minimized.

The sources of radiated EMI can be broadly classified into three types:

1. Hand-held portable transceivers (transmitter-receivers) with the antenna mounted directly on the transmitting unit. Examples include: citizens band (CB) radios, walkie talkies, security, fire, and police transceivers, cellular telephones and other personal communication devices.

NOTE: Some cellular telephones and similar transmit signal while they are ON, even when not being used.

2. Medium-range mobile transceivers, such as those used in police cars, fire trucks, ambulances and taxis. These usually have the antenna mounted on the outside of the vehicle.

4. ELECTROMAGNETIC INTERFERENCE

3. Long-range transmitters and transceivers, such as commercial broadcast transmitter (radio and TV broadcast antenna towers) and amateur (HAM) radios.

Note: Other types of hand-held devices, such as cordless phones, laptop computers, AM/FM radios, TV sets, CD players, and cassette players, and small appliances, such as electric shavers and hair dryers, so far as we know, are not likely to cause EMI problems to your powered vehicle.

➤ **POWERED VEHICLE ELECTROMAGNETIC INTERFERENCE (EMI)**

Because EM energy rapidly becomes more intense as one moves closer to the transmitting antenna (source), the EM fields from hand-held radio wave sources (transceivers) are of special concern. It is possible to unintentionally bring high levels of EM energy very close to the powered vehicle's control system while using these devices. This can affect powered vehicle movement and braking. Therefore, the warnings listed below are recommended to prevent possible interference with the control system of the powered vehicle.

➤ **WARNINGS**



Electromagnetic interference (EMI) from sources such as radio and TV stations, amateur radio (HAM) transmitters, two-way radios, and cellular phones can affect electric powered vehicles and motorized vehicles. Following the warnings listed below should reduce the chance of unintended brake release or electric powered vehicle movement which could result in serious injuries.

1. Do not operate hand-held transceivers-receivers, such as citizens band (CB) radios, or turn ON personal communication devices, such as cellular phones, while the electric powered vehicle is turned ON.
2. Be aware of nearby transmitters, such as radio or TV stations, and try to avoid coming close to them.
3. If unintended movement or brake release occurs, turn the powered vehicle OFF as soon as it is safe.
4. Be aware that adding accessories or components, or modifying the electric powered vehicle, may make it more susceptible to EMI (Note: There is no easy way to evaluate their effect on the overall immunity of the powered vehicle)
5. Report all incidents of unintended movement or brake release to the powered vehicle manufacturer, and note whether there is a source of EMI nearby

5. PARTS INTRODUCTION

➤ PARTS DESCRIPTION, DK S600

1.Speed Dial	2. Rear Mirror	3. Control Panel
4.Tiller	5. Basket	6. Head Lamp
7. Signal Lamp	8. Front bumper	9. Head – Rest
10. Arm Rest	11. Seat	12. Seat Angle Adjust Bar
13. Cup/Bottle Holder	14. Rear Light Package	15. On-board charger
16. Freewheel Mode Lever	17. Anti-tip wheels	



6. OPERATION

6.1 OPERATION PANEL

1. Speed Dial
2. Horn Button
3. Signal Light Switch
4. Key Hole/Head Light / Rear Light Switch
5. Speed Control Lever
6. Battery Condition Meter
7. Warning light switch
8. High Speed / Low speed Switch



6.2 HOW TO OPERATE YOUR VEHICLE

➤ Power switch

Turn the key to power ON or OFF

OFF: Power is turned off

ON: Power is turned on without light on

LIGHT: Both power and light are turned on



➤ Speed Dial

Allows the user to set the speed to a comfortable and safe level.



➤ Forward, Reverse, and Braking

- ◆ Pull the speed control lever down with your right hand and the vehicle will move forward.
- ◆ Pull the speed control lever down with your left hand and the vehicle will move backward.
- ◆ Release the speed control lever freely while in either forward or reverse direction, the electromagnetic brake in the motor will be activated and the vehicle will stop

6. OPERATION

➤ Signal Lamp Switch

- ◆ The front and rear lights will flash in the direction of the turn signal indicated by the user.
- ← left signal lamp
- Right signal lamp
- ◆ Push the turn signal again and it will turn the signal off.



➤ Horn Button

The warning horn will sound when pressed.

■ Key Hole/ Power Switch/ Light Switch

Turn the key to “ON” position to start the vehicle. Turn key to “LIGHT” will turn on both the front and rear lights.





WARNING

When going down a hill, NEVER set to freewheel mode. The electromagnetic brake will not activate

➤ Seat

- The seat can be turned a full 360 degrees in 90 degree increments using the rotation lever. (1)
- Seat height can be adjusted by removing the adjustment pin (3) and loosening the adjustment knob (4). To access the pin and the adjustment knob you must first remove the battery cover (2) using a Phillips screwdriver.
- The seat can slide forwards and backwards. Locate the adjustment lever on the side of the seat. (5) Lift to adjust.
- The seat can be put in a reclined position by using the lever located on the side of the seat. (6)





6. OPERATION

➤ Battery Condition Meter

- ◆ When the power switch is turned on, the battery condition meter will light up with 5 LED lights indicating the remaining power of battery.
All 5 LED lights: full power. 1 LED light with beeper warning: out of power.
- ◆ The remaining power indicated by battery meter will vary by the driving time and the way the vehicle is driven. Repeated starting, stopping and climbing will consume the power more quickly.



i SUGGESTION

1. You are strongly recommended to charge the battery immediately when the battery meter shows less than 2 LED lights.
2. After charging or replacing the batteries, drive the vehicle for 2-3 minutes to make sure the batteries hold a charge.
3. In winter, the batteries may respond slower and the traveling distance may be less than in the summer.
4. The battery indicator, while driving up an incline, may move up and down. This is normal.
5. Under normal circumstances, it is natural for the capacity of the batteries to decay as time passes which can result in a short travel distance compared to a brand new battery. When you find the travel distance is about only 50% of a brand new battery's traveling distance, it is time for replacement batteries. Contact your dealer and ask for replacement batteries. If you continue to use the old battery less travel distance will be expected.
6. Travel distances will be shortened when driving frequently on a slope or ground with big dips, due to more power consumption.


i SUGGESTION

When driving indoors or in a busy area, for your own safety, Be sure to set the switch to “L” which has lower max speed. You are recommended to set the switch to “H” only at a wide open area.

6. OPERATION

6.3 HOW TO SET TO FREEWHEEL MODE

- ◆ Drive on mode: Push lever down completely, the vehicle can be driven by motor power.
- ◆ Freewheel Mode: Lift the lever up completely, the vehicle can now be moved manually.



➤ Tiller Adjustment

The Tiller can be adjusted between 0° ~ 80° by following steps.

1. Pull the adjusting lever on the tiller.
2. Adjust the tiller to the desired angle.
3. Let the lever go, and tiller will automatically fix in the desired angle.



WARNING

- ◆ Do not set the vehicle in freewheel mode when driving on an incline or decline. This could lead to the user being unable to control the speed of the vehicle.
- ◆ To protect yourself, while you are driving down a steep decline (over 13°), the vehicle will not be driven automatically and electromagnetic brake will activate. Turn the power off and on again to restart your vehicle.

6. OPERATION

➤ Starting and Driving

1. Make sure the seat is installed properly.
2. Make sure the tiller has been secured properly.
3. Make sure the armrest has been returned to the lowest position.
4. Turn the power switch to “ON.” If necessary, turn on the headlights.
5. Check battery indicator to see if there is enough power for your trip. If you have any doubts, please charge the battery before leaving.
6. Set the proper speed.
7. Be sure the speed lever is in working condition.
8. Make sure the electromagnetic brake is in working condition.

! ATTENTION

1. Do not push both Right Hand & Left Hand side speed control levers simultaneously. This might lead you to be unable to control your vehicle.
2. Do not turn the power switch to OFF while driving, as this will lead to an emergency stop and possible injury.
3. Do not set to high speeds while driving indoors.
4. Do not adjust speed dial while driving, the sudden change in speed may endanger you and cause the failure of your vehicle.
5. Do not place magnetic devices near the area of the operation handle as this could affect the safe operation of your vehicle.
6. Be careful while driving in heavy traffic or crowded areas.
7. While in reverse, be aware of people or objects around you.

Speed Scale	Circumstances of Recommend speed
Left side	<ul style="list-style-type: none">◆ Speed equivalent to walking slowly◆ Use while on an incline or decline◆ Use while driving indoors or in a narrow space
Between left and right	<ul style="list-style-type: none">◆ Speed equivalent to walking normally◆ Use while driving outdoors on flat ground
Right side	<ul style="list-style-type: none">◆ Speed equivalent to walking quickly◆ Use while driving in an open space.

6. OPERATION

➤ Stopping

1. Release the speed control lever completely, and the vehicle will activate the electromagnetic brake and stop.
2. Turn the power switch to “OFF.” Then pull out the key.

! ATTENTION

- ◆ The stopping distance will vary with your speed.
- ◆ While parking your vehicle, be sure to park on flat ground and then turn the power to “OFF” before you get off.

7. CHARGING AND BATTERY

7.1 CHARGING THE BATTERY

Be sure to follow all procedures below

- Turn the power switch to “OFF”.
- Open the charging socket cap in the rear lighting package
- Connect the charger’s cord onto the power outlet. The length of cord is limited so DO NOT pull out the cord beyond the yellow mark otherwise the cord will not be able to roll back itself.
- Charger LED light will be on at the beginning of charging. The charging duration is 8 hours up to 10 hours, based on the status of the battery and ambient temperature.
- The red LED will change to green when the charging is complete.
- Disconnect the cord; pull out the cord slightly and let the cord roll back.



SUGGESTION

1. Do not disconnect the charger cord if the charging is not complete. The battery life will be seriously shortened or decayed if the battery is repeatedly used without being fully charged.
2. The vehicle will indicate a full charge by the charging light turning from red to green.
3. If the vehicle is not plugged in, the batteries will slowly discharge over time.
4. If you do not use your vehicle for long periods of time, it should be charged at least every three months to keep the battery in a full charge status.

ATTENTION: Please follow the rules below to avoid accidents while charging.

1. Never disassemble or modify the plug.
2. Please charge in a well ventilated space where is not directly exposed to the sunlight. Do not charge in an area which it is humid, under rainfall, or morning dew.
3. Do not cover the vehicle with any waterproof cloth or other objects while charging.
4. Do not charge your vehicle in temperatures less than 14 degrees F or higher than 122 degrees F, as the charger may not work well and the battery may be damaged.

7. CHARGING AND BATTERY

WARNING

1. Keep away from flammable objects while charging, as it may lead to fire or explosion of battery.
2. As the battery may produce hydrogen, do not smoke while charging. Always charge your batteries in a well-ventilated space.
3. Never connect or disconnect the charging cord with wet hands while charging. Do not connect or disconnect the charging cord when it is wet, it may lead to an electric shock.

7.2 CHARGING CORD AND PLUG

- ◆ Never pull out the charging cord beyond the **yellow mark**. Otherwise the cord won't automatically roll back. Please note pulling out beyond the **red mark** will damage the cord and charger.



WARNING

Park the vehicle near the electric outlet and slowly pull out the charging cord. If the charging cord is broken, please **DO NOT** charge the vehicle until a new charging cord is obtained.

7.3 BATTERY

- ◆ Do not expose the battery to temperatures below 14 degrees F or above 122 degrees F when charging or storing the vehicle. Exceeding the specified temperature ranges can lead to freezing or over heating of the batteries. Such case will damage the batteries and shorten their lifetime.
- ◆ The battery used on your vehicle is sealed lead acid battery, maintenance-free, so it is not necessary to change the battery liquid or refill the liquid.

WARNING

Do not open the battery seal cap at any time.

7. CHARGING AND BATTERY

➤ **Cleaning the battery**

If the battery is contaminated by dirty water, battery acid, dust, or other factors, the battery will discharge quickly. Follow the following steps to clean the battery.

1. Turn the vehicle power switch to OFF.
2. Remove the seat and dust cover.
3. Remove the shroud and unplug the terminal of the tail light and signal lights.
4. Use a clean cloth to wipe off the soiled area.
5. Take out the battery.
6. Clean the battery with a clean cloth. If the terminal is covered with white powder, please clean it with warm water.



SUGGESTION

If necessary, contact your B Mobile service center for cleansing and replacing of the battery.



SUGGESTION

1. Make sure the terminals are installed properly then put the cover back on
2. Do not use the battery to charge telecom equipment or other items.
3. Always replace both batteries simultaneously.

8. INSPECTION AND MAINTENANCE

8.1 DAILY CHECKING

Check following items before driving. If you find any thing abnormal, go to your B Mobile service center for further inspection before using it.

Item	Inspection Content
Handle bar	<ul style="list-style-type: none">◆ Are they loose?◆ Can they be turned left, right smoothly?
Speed Control Dial	<ul style="list-style-type: none">◆ Can it be adjusted freely and function well?
Speed Control Lever	<ul style="list-style-type: none">◆ Does the vehicle will move when the lever is pressed down?◆ Does the vehicle stop when the lever is released completely?
Motor	<ul style="list-style-type: none">◆ Is there any abnormal noise from the motor?◆ Does the electromagnetic brake work properly?
Freewheel Mode	<ul style="list-style-type: none">◆ Does the freewheel mode lever work properly?
Battery Condition meter	<ul style="list-style-type: none">◆ Is the light on when power is switched on?◆ Is the remaining power enough for your trip?
Horn	<ul style="list-style-type: none">◆ Does the horn work well?
Turn signal	<ul style="list-style-type: none">◆ Does the turn signal work?
Headlight	<ul style="list-style-type: none">◆ Does the headlight work?
Taillight	<ul style="list-style-type: none">◆ Does the taillight work?
Seat	<ul style="list-style-type: none">◆ Can the seat be turned and adjusted smoothly?
Rearview Mirror	<ul style="list-style-type: none">◆ Are they adjusted properly or broken?
Tires	<ul style="list-style-type: none">◆ Are there any cracks or other damage to the tires?◆ Are the tires properly inflated?◆ Check the tires tread depth.
Other	<ul style="list-style-type: none">◆ Are there any abnormal noises?◆ Is there an oil leakage from transmission box?

! ATTENTION

Go to Bladez Health Products Corporation service center for inspection and maintenance if you find anything wrong.

8. INSPECTION AND MAINTENANCE

8.2 REGULAR CHECKING RECORD

To make sure your vehicle remains in good condition, go to your B Mobile service center regularly for maintenance and record it accordingly every six months after purchasing. (Fee Required)

SUGGESTION

Even if you don't use the vehicle for a long time, the vehicle should still be maintained regularly.

I – Inspection A-Adjustment R-Repair RP- Replacement T - Tightening
C – Clean L - Add Oil

Item \ Period		1 mo	6 mo	12 mo	18 mo	24 mo	30 mo
Operate Panel	Forward/Reverse of speed lever						
	The function of switches						
	Connection of Socket / Plug is good?						
Motor	Operation and Noise						
	Function of Electromagnetic Brakes						
	Connection of Socket / Plug is good?						
Transmi -ssion	Freewheel mode lever						
	Any oil leakage from transmission box						
Battery	Looseness or corrosion on terminals						
	Connection of Socket / Plug is good						
	Any broken on electrical cord						
Charger	Charging function and LED light						
	Any broken on Socket / Plug / Cord						
Seat	Angle adjustment bar						
	Seat post looseness						
Tiller	Right /Left Turn of tiller						
	Adjustment of Tiller angle						
Tire	Crack or damage of tire surface						
	Tread depth enough or not						
	Abnormal wearing						
	Tire Pressure						

8. INSPECTION AND MAINTENANCE

Item \ Period		1 mo	6 mo	12 mo	18 mo	24 mo	30 mo
Wheel Rim	Tightness of bolts and nuts						
	Wear or deformation around tires						
Shock Absorber	Tightness of bolts and nuts						
	Damage to shock absorbers						
Inspector	Checked By						
	Date Checked						

8.3 FUSE AND TIRES

➤ Fuse

There are two fuses on your vehicle. One is located inside the handle bar cover; the other is near the battery. If the power switch is “ON” and the battery condition meter is not lit, check the fuse.

SUGGESTION

Contact your Bladez Health Product Corporation service center for inspecting or replacing the fuse, since the tiller shroud has to be removed first to replace the fuse.

➤ Tires

The condition of the tires depends on how you drive and use your vehicle.

- ◆ When the tread depth is below 0.5mm, it can easily lead to wheel slippage, making braking distance longer. Replace the tires as soon as possible.
- ◆ The tire pressure shall be kept at about 60 PSI (approx 4.2KG/cm) for the best travel range and usage.

8. INSPECTION AND MAINTENANCE

8.4 MAINTENANCE

- ◆ You must maintain the vehicle frequently if you constantly drive on grass, sand, or gravel road.
- ◆ Do not use water, oil or other chemical solution to clean your vehicle. DO NOT spray the vehicle with water in order to prevent damaging the electronic components and PC circuit board. Please clean the vehicle by wiping it with either dry or moist cloth.
- ◆ Please take the vehicle to authorized B Mobile service center for repairs and adjustments. Improper adjustments could lead to accidents and vehicle malfunction.
- ◆ Please use soft and dry cloths to keep your vehicle tidy. Please use moderate or mild detergent to clean the vehicle.

! ATTENTION

When conducting the maintenance of your vehicle, please turn the power switch OFF and remove the charger cords.

i SUGGESTION

- ◆ Do not use gasoline, solvents or vaporizing solution, the shroud may be deformed or damaged.
- ◆ Do not use wax.

8.5 STORAGE

Make sure the vehicle is stored under the following circumstances:

- ◆ Make sure the Power switch is turned to “OFF.”

9. TROUBLESHOOTING

➤ **Vehicle Turns On But Doesn't Move**

- ◆ May still be in freewheel mode.
- ◆ Make sure the lever is completely pushed down to disengaged freewheel mode, then turn off the vehicle and turn it back on.
- ◆ This problem is indicated by the light blinking five times.

➤ **Vehicle Moves Very Slow**

- ◆ You might be exceeding the maximum weight limit of the vehicle, refer to “Maximum Weight Limit” and try again.
- ◆ Your vehicle Batteries might be low. Check the battery meter on the front display and charge if necessary.
- ◆ Also check for obstructions on or around the vehicle of in your pathway.
- ◆ Check high/low and adjustable speed knob.

➤ **Vehicle Seems To Be Getting No Power**

- ◆ Check the fuses which are located in the handlebar cover and near the batteries.
- ◆ Check the Batteries connections by removing the rear shroud, also inspect the wires to see if they have been pinched or frayed.
- ◆ If the vehicle hasn't been charged regularly, the batteries may be completely discharged, Charge the vehicle and try again.

9. TROUBLESHOOTING

DC Status Indicator Flash Code	Fault Description	Impact on Controller	Notes
1	Battery needs recharging	Reduced drive is available	The battery voltage has dropped below 23.3 volts in neutral. Recharge batteries.
2	Battery voltage too low	Drive Inhibited	The battery charge has dropped past 16.5 volts and is not sufficient to allow safe driving.
3	Battery voltage too high	Drive inhibited. If the vehicle is traveling, a soft stop occurs.	The battery charge has exceeded 32 volts. Check the battery conditions and the connections.
4	Current limit time out	Drive inhibited	The motor current has reached too high a value. Check the condition of the motor and loom(s). Contact your service agent
5	Freewheel engaged	Drive inhibited	Check brake connections. Check the condition of the motor and the loom(s). Contact your service agent.
6	Speed control out of neutral	Drive inhibited	Return speed to neutral and reset system. Readjust the speed lever neutral if necessary.
7	Speed control fault	Drive inhibited	Check speed wiring for open or short circuits. Check speed set-up. Contact your service agent.
8	Motor fault	Drive inhibited	Contact your service agent
9	Internal controller fault	Drive inhibited	Contact your service agent

10. SPECIFICATIONS

	DK S600
Model Name	B*Mobile – President
Packing Size	64.2"x28.0"x26.3" (1630mmx710mmx670mm)
Overall Height (without seat)	38" (965.20 mm)
Overall Length	47" (1193.8 mm)
Overall Width	26.5" (673.10 mm)
Motor Output	700W
Battery	12Vx50AHx2PCs
Charger	24V 3AMP (Automatic Type) on-board
Front Wheel	3.3"x10.2" (85mmx260mm) pneumatic tire x 2 PCs
Rear Wheel	8"x 13.2" (200mmx335mm) pneumatic tire x 2 PCs
Break System	Intelligent regenerative electromagnetic brake
Net Weight (with battery)	206.8 lbs (93.8 kg)
Net Weight (without battery)	135 lbs (61.4 kg)
Maximum Loading	450 lbs (204.1kg)
Slope Grade Ability	13 degrees
Continuous Traveling Distance	28.0 miles (45 km)
Maximum Speed	7.5 mph
Controller	Dynamic Rhino 110 AH
Turning Radius	42.3" (1074.4mm)
Maximum Ground Clearance	2.8" (70mm)
Front Suspension	YES
Rear Suspension	YES

Remark -- The manufacturer reserves the right to modify the specifications if needed. The final specifications are subject to any change when you purchase the vehicle at the shop.

Maximum driving distance is based on the conditions of an ambient temperature, a 150 lbs driver and brand-new fully charged battery by a constant driving speed at 3.7 mph with 70% of battery power being discharged.

11. WARRANTY

11.1 VSN (VEHICLE SERIAL NUMBER)

To ensure after-sale service and warranty service, please fill in the vehicle serial number that is stuck on the back right hand side of the frame.

Model	DK S600
VSN	

11.2 WARRANTY POLICY

11.2.1 Contents and Duration of Warranty

➤ Contents of Warranty

1. For one year from the date of purchase, Bladez Health Products Corp will replace at our option free of charge, any part or electronic component found to be defective due workmanship and/or material upon the examination by a certified Bladez Health Product Corp. service center.
2. Warranty service can be performed by authorized B Mobile service centers only. Do not return faulty parts to Bladez Health Products Corp. with out prior consent.
3. Transportation, repair, or replacement costs will be the purchaser's responsibility.
4. Replaced parts will be the properties of Bladez Health Products Corp. Do not return faulty part to Bladez Health Products Corporation without prior consent.

NOTE: Some states do not allow limitations on how long an implied warranty lasts or does not allow the exclusion of limitation of incidental or consequential damages.

➤ Warranty duration and parts covered by the warranty

Parts and labor	1 year from the date of purchase
General Items	1 year from the date of purchase
Battery	Warranted by the battery manufacturer. The battery is not warranted by Bladez Health Products Corp.

11. WARRANTY

11.2.2 Items Not Covered By Warranty

➤ Items listed below

1. Consumable parts such as brake shoe or pad, clutch shoe and lining, bulbs, fuses, upholstery and seating, brake cable, oil seal packing and gasket, screw and washer, lubricant oil and grease, carbon brush inside motor, battery liquid, and footrest carpet etc.
2. A sensation which does not affect the quality and function of the vehicle, such as noise or vibration.
3. Wear and tear resulting from usage as well as circumstances resulting from the lapse of time. (Natural aging and deterioration of paint, electroplating, plastic parts, etc)
4. Damage due to lack of maintenance or use of improper operation or storage.
5. Fees for procedural inspections, adjustments, adding oil, cleaning, and other maintenance.
6. Fees for the regular inspections and maintenance.
7. Any fees incurred as a result of the warranty repair, losses or compensation because of the inability to use the electric vehicle (telephone use, shipping, car rental, travel fees, etc)

➤ Problems caused by the following reasons are not covered:

1. Failure to carry out the regular inspection and maintenance recommended by this owner's manual.
2. Improper storage.
3. Tampering by changing or adding components to the vehicle which is not approved by Bladez Health Products Corp.
4. Damage caused by not using Bladez Health Products Corp. genuine parts or any parts not recommended by Bladez Health Products Corp.
5. Abuse or misuse of the vehicle such as driving over curbs or overloading with passengers or goods.
6. Driving in places that the vehicle should not normally be driven or used.
7. Damage or surface corrosion from environment such as airborne fallout (chemicals, tree sap, etc.) stones, hail, windstorm, lightning, flood, etc.

11.2.3 The responsibility of user

1. Proper use of the vehicle in accordance with this Owner's Manual.
2. The carry out of daily inspections
3. The carry out of the regular maintenance set by this company.
4. Keep a record of vehicle inspection and maintenance in this Owner's Manual.

11. WARRANTY

11.2.4 Warranty effectiveness

The Product registration Card must be completely filled out and stamped by the Bladez Health Products Corporation Dealer where it was purchased.

11.2.5 Warranty rights are non-Transferable

The Bladez Health Products Corp warranty is non-Transferable. The warranty will only be honored to the original purchaser.

! ATTENTION: Implied warranties, including those of merchantability and fitness for a particular purpose, are limited to one year from the date of original purchase and to the extent permitted by law. Any and all implied warranties are excluded. This is the exclusion remedy. Liabilities for consequential damages under any and all warranties are excluded.